

产品名称: Nkx-2.8 Rabbit Polyclonal Antibody
产品货号: APRab14732



产品概述 (Summary)

产品名称 (Production Name)	Nkx-2.8 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NKX2-8
别名 (Alternative Names)	
基因 ID (Gene ID)	26257.0
蛋白 ID (SwissProt ID)	Q8IUT7.The antiserum was produced against synthesized peptide derived from human NKX28. AA range:20-69

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	26kDa

研究背景 (Background)


产品名称: Nkx-2.8 Rabbit Polyclonal Antibody
产品货号: APRab14732



The protein encoded by this gene is a homeobox-containing developmental regulator associated with liver development. The encoded protein binds to the alpha-fetoprotein (AFP) gene promoter and increases the expression of AFP. This gene is overexpressed in some lung cancers and is linked to poor patient survival, possibly due to its resistance to cisplatin. This gene is aberrantly methylated in pancreatic cancer, deleted in squamous cell lung carcinomas, and acts as a tumor suppressor in esophageal cancer. Mutations in this gene may also be a cause of neural tube defects. [provided by RefSeq, Dec 2015],similarity:Belongs to the NK-2 homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of NKX28 Antibody. The lane on the right is blocked with the NKX28 peptide.

注意事项 (Note)

For research use only .