

产品名称: Nkx-6.1 Rabbit Polyclonal Antibody
产品货号: APRab14737



产品概述 (Summary)

产品名称 (Production Name)	Nkx-6.1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NKX6-1
别名 (Alternative Names)	NKX6-1; NKX6A; Homeobox protein Nkx-6.1; Homeobox protein NK-6 homolog A
基因 ID (Gene ID)	4825.0
蛋白 ID (SwissProt ID)	P78426.Synthesized peptide derived from Nkx-6.1 . at AA range: 180-260

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

产品名称: Nkx-6.1 Rabbit Polyclonal Antibody
产品货号: APRab14737




In the pancreas, NKX6.1 is required for the development of beta cells and is a potent bifunctional transcription regulator that binds to AT-rich sequences within the promoter region of target genes lype et al. (2004) [PubMed 15056733]. [supplied by OMIM, Mar 2008], function: May be important for control of islet development and/or regulation of insulin biosynthesis., similarity: Contains 1 homeobox DNA-binding domain., tissue specificity: Pancreatic beta cells.,

研究领域 (Research Area)

Maturity onset diabetes of the young;

图片 (Image Data)



Western blot analysis of NKX61 Antibody. The lane on the right is blocked with the NKX61 peptide.

注意事项 (Note)

For research use only .