

产品名称: NOC2L Rabbit Polyclonal Antibody
产品货号: APRab14774



产品概述 (Summary)

产品名称 (Production Name)	NOC2L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOC2L
别名 (Alternative Names)	NOC2L; NIR; Nucleolar complex protein 2 homolog; Protein NOC2 homolog; NOC2-like protein; Novel INHAT repressor
基因 ID (Gene ID)	26155.0
蛋白 ID (SwissProt ID)	Q9Y3T9.The antiserum was produced against synthesized peptide derived from human NOC2L. AA range:602-651

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	85kDa

产品名称: NOC2L Rabbit Polyclonal Antibody
产品货号: APRab14774



研究背景 (Background)


NOC2 like nucleolar associated transcriptional repressor(NOC2L) Homo sapiens Histone modification by histone acetyltransferases (HAT) and histone deacetylases (HDAC) can control major aspects of transcriptional regulation. NOC2L represents a novel HDAC-independent inhibitor of histone acetyltransferase (INHAT) (Hublitz et al., 2005 [PubMed 16322561]).[supplied by OMIM, Mar 2008],similarity:Belongs to the NOC2 family.,

研究领域 (Research Area)


Transcription; Co-factors; Epigenetics and Nuclear Signaling; Chromatin Modifying Enzymes; Acetylation; HAT; Cell cycle; Cell Cycle Inhibitors; p53;

Cell Biology; Apoptosis; Intracellular; p53 Pathway

图片 (Image Data)



Western blot analysis of lysates from MCF-7 cells, using NOC2L Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody was pre-absorbed by immunogen peptide.

产品名称: NOC2L Rabbit Polyclonal Antibody
产品货号: APRab14774



注意事项 (Note)

For research use only .