

**产品名称: Nodal Rabbit Polyclonal Antibody**  
**产品货号: APRab14780**



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | Nodal Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA              |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                  |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                   |
|------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | NODAL                                                                                                                             |
| 别名 (Alternative Names) | NODAL; Nodal homolog                                                                                                              |
| 基因 ID (Gene ID)        | 4838.0                                                                                                                            |
| 蛋白 ID (SwissProt ID)   | Q96S42.The antiserum was produced against synthesized peptide derived from the C-terminal region of human NODAL. AA range:294-343 |

## 产品应用 (Application)

|                          |                                                                        |
|--------------------------|------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 40kDa                                                                  |

## 研究背景 (Background)

**产品名称: Nodal Rabbit Polyclonal Antibody**  
**产品货号: APRab14780**




This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate the mature protein, which regulates early embryonic development. This protein is required for maintenance of human embryonic stem cell pluripotency and may play a role in human placental development. Mutations in this gene are associated with heterotaxy, a condition characterized by random orientation of visceral organs with respect to the left-right axis. [provided by RefSeq, Aug 2016],disease:Defects in NODAL may be a cause of situs ambiguus [MIM:601265].,function:Essential for mesoderm formation and axial patterning during embryonic development.,similarity:Belongs to the TGF-beta family.,subunit:Homodimer; disulfide-linked.,


## 研究领域 (Research Area)

TGF-beta;

## 图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100



Western blot analysis of HeLa lysis using NODAL antibody. Secondary antibody was diluted at 1:20000

## 注意事项 (Note)

For research use only .

产品名称: Nodal Rabbit Polyclonal Antibody  
产品货号: APRab14780

---

