

产品名称: NOL7 Rabbit Polyclonal Antibody
产品货号: APRab14790



产品概述 (Summary)

产品名称 (Production Name)	NOL7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOL7 C6orf90 NOP27
别名 (Alternative Names)	
基因 ID (Gene ID)	51406.0
蛋白 ID (SwissProt ID)	Q9UMY1.Synthesized peptide derived from human protein . at AA range: 190-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	28kDa

研究背景 (Background)

The protein encoded by this gene localizes to the nucleolus, where it maintains nucleolar structure and cell growth rates.


产品名称: NOL7 Rabbit Polyclonal Antibody
产品货号: APRab14790



The encoded protein also functions as a tumor suppressor and regulator of angiogenesis. The RB tumor suppressor gene recruits transcription factors to this gene and positively regulates its expression. Two transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Dec 2015],sequence caution:Contaminating sequence. Potential poly-A sequence.,tissue specificity:Expressed in numerous tissues. Particularly prevalent in the adrenal gland, thyroid gland, heart and skeletal muscle.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .