

产品名称: NOM1 Rabbit Polyclonal Antibody
产品货号: APRab14793



产品概述 (Summary)

产品名称 (Production Name)	NOM1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOM1
别名 (Alternative Names)	NOM1; C7orf3; Nucleolar MIF4G domain-containing protein 1; SGD1 homolog
基因 ID (Gene ID)	64434.0
蛋白 ID (SwissProt ID)	Q5C9Z4. The antiserum was produced against synthesized peptide derived from human NOM1. AA range: 661-710

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: NOM1 Rabbit Polyclonal Antibody
产品货号: APRab14793



Proteins that contain MIF4G (middle of eIF4G (MIM 600495)) and/or MA3 domains, such as NOM1, function in protein translation. These domains include binding sites for members of the EIF4A family of ATP-dependent DEAD box RNA helicases (see EIF4A1; MIM 602641) (Simmons et al., 2005 [PubMed 15715967]).[supplied by OMIM, Mar 2008],similarity:Belongs to the CWC22 family.,similarity:Contains 1 MI domain.,similarity:Contains 1 MIF4G domain.,subunit:May interact with EIF4A1, EIF4A2 and EIF4A3.,tissue specificity:Expressed in heart and skeletal muscle.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using NOM1 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .