

产品名称: Nov Rabbit Polyclonal Antibody
产品货号: AP Rab14809



产品概述 (Summary)

产品名称 (Production Name)	Nov Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NOV NOV; CCN3; IGFBP9; NOVH; Protein NOV homolog; NovH; CCN family member 3; Insulin-like growth factor-binding protein 9; IBP-9; IGF-binding protein 9; IGFBP-9; Nephroblastoma-overexpressed gene protein homolog
别名 (Alternative Names)	
基因 ID (Gene ID)	4856.0
蛋白 ID (SwissProt ID)	P48745.The antiserum was produced against synthesized peptide derived from the Internal region of human NOV. AA range:221-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	39kDa

产品名称: Nov Rabbit Polyclonal Antibody
产品货号: APRab14809




研究背景 (Background)


The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CCN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development. [provided by RefSeq, Feb 2009],function:Immediate-early protein likely to play a role in cell growth regulation.,similarity:Belongs to the CCN family.,similarity:Contains 1 CTCK (C-terminal cystine knot-like) domain.,similarity:Contains 1 IGFBP N-terminal domain.,similarity:Contains 1 TSP type-1 domain.,similarity:Contains 1 VWFC domain.,subunit:Interacts with FBLN1.,tissue specificity:Increased expression in Wilms tumor of the stromal type.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysate from MCF-7 cells, using NOV Antibody.



Western Blot analysis of MCF-7 cells using Nov Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

产品名称: Nov Rabbit Polyclonal Antibody
产品货号: APRab14809



Immunohistochemical analysis of paraffin-embedded human-skin, antibody was diluted at 1:100

注意事项 (Note)

For research use only .