

产品名称: NPDC-1 Rabbit Polyclonal Antibody
产品货号: APRab14834



产品概述 (Summary)

产品名称 (Production Name)	NPDC-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPDC1
别名 (Alternative Names)	NPDC1; Neural proliferation differentiation and control protein 1; NPDC-1
基因 ID (Gene ID)	56654.0
蛋白 ID (SwissProt ID)	Q9NQX5.The antiserum was produced against synthesized peptide derived from human NPDC1. AA range:189-238

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)


产品名称: NPDC-1 Rabbit Polyclonal Antibody
产品货号: APRab14834




function:Suppresses oncogenic transformation in neural and non-neural cells and down-regulates neural cell proliferation. Might be involved in transcriptional regulation.,similarity:Belongs to the NPDC1/cab-1 family.,tissue specificity:Strongly expressed in adult brain; especially in hippocampus, frontal lobe and temporal lobe.,function:Suppresses oncogenic transformation in neural and non-neural cells and down-regulates neural cell proliferation. Might be involved in transcriptional regulation.,similarity:Belongs to the NPDC1/cab-1 family.,tissue specificity:Strongly expressed in adult brain; especially in hippocampus, frontal lobe and temporal lobe.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HepG2 and LOVO cells, using NPDC1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using NPDC1 antibody.

产品名称: NPDC-1 Rabbit Polyclonal Antibody
产品货号: APRab14834



Western Blot analysis of various cells using NPDC-1 Polyclonal Antibody

注意事项 (Note)

For research use only .