

产品名称: NPY1-R Rabbit Polyclonal Antibody
产品货号: AP Rab14859



产品概述 (Summary)

产品名称 (Production Name)	NPY1-R Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NPY1R
别名 (Alternative Names)	NPY1R; NPYR; NPY1; Neuropeptide Y receptor type 1; NPY1-R
基因 ID (Gene ID)	4886.0
蛋白 ID (SwissProt ID)	P25929. The antiserum was produced against synthesized peptide derived from human NPY1R. AA range: 221-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	48kDa

研究背景 (Background)

产品名称: NPY1-R Rabbit Polyclonal Antibody
产品货号: AP Rab14859




This gene belongs to the G-protein-coupled receptor superfamily. The encoded transmembrane protein mediates the function of neuropeptide Y (NPY), a neurotransmitter, and peptide YY (PYY), a gastrointestinal hormone. The encoded receptor undergoes fast agonist-induced internalization through clathrin-coated pits and is subsequently recycled back to the cell membrane. Activation of Y1 receptors may result in mobilization of intracellular calcium and inhibition of adenylate cyclase activity. [provided by RefSeq, Aug 2013],function:Receptor for neuropeptide Y and peptide YY. The rank order of affinity of this receptor for pancreatic polypeptides is NPY > [Pro-34] PYY, PYY and [Leu-31, Pro-34] NPY > NPY (2-36) > [Ile-31, Gln-34] PP and PYY (3-36) > PP > NPY free acid.,similarity:Belongs to the G-protein coupled receptor 1 family,


研究领域 (Research Area)

Neuroactive ligand-receptor interaction;

图片 (Image Data)




Western blot analysis of lysates from NIH/3T3 cells, using NPY1R Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from 293 cells using NPY1R antibody.

产品名称: NPY1-R Rabbit Polyclonal Antibody
产品货号: AP Rab14859



Western Blot analysis of 3T3 cells using NPY1-R Polyclonal Antibody

注意事项 (Note)

For research use only .