

产品名称: NRBF-2 Rabbit Polyclonal Antibody
产品货号: APRab14885



产品概述 (Summary)

产品名称 (Production Name)	NRBF-2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NRBF2
别名 (Alternative Names)	NRBF2; COPR; Nuclear receptor-binding factor 2; NRBF-2; Comodulator of PPAR and RXR
基因 ID (Gene ID)	29982.0
蛋白 ID (SwissProt ID)	Q96F24.The antiserum was produced against synthesized peptide derived from human NRBF2. AA range:140-189

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

产品名称: NRBF-2 Rabbit Polyclonal Antibody
产品货号: APRab14885




研究背景 (Background)


function: May modulate transcriptional activation by target nuclear receptors. Can act as transcriptional activator (in vitro)., subunit: Interacts with PPARA, PPARD and PPARG. Interacts with RARA, RARG and RXRA in the presence of bound ligand., tissue specificity: Detected in keratinocytes, liver and placenta., function: May modulate transcriptional activation by target nuclear receptors. Can act as transcriptional activator (in vitro)., subunit: Interacts with PPARA, PPARD and PPARG. Interacts with RARA, RARG and RXRA in the presence of bound ligand., tissue specificity: Detected in keratinocytes, liver and placenta.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from A549 cells, using NRBF2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7 cells using NRBF2 antibody.

产品名称: NRBF-2 Rabbit Polyclonal Antibody
产品货号: APRab14885



Western Blot analysis of various cells using NRBF-2 Polyclonal Antibody

注意事项 (Note)

For research use only .