

产品名称: NRIP3 Rabbit Polyclonal Antibody
产品货号: APRab14898



产品概述 (Summary)

产品名称 (Production Name)	NRIP3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NRIP3
别名 (Alternative Names)	NRIP3; C11orf14; Nuclear receptor-interacting protein 3; Sarcoma antigen NY-SAR-105
基因 ID (Gene ID)	56675.0
蛋白 ID (SwissProt ID)	Q9NQ35.The antiserum was produced against synthesized peptide derived from human NRIP3. AA range:192-241

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	27kDa

产品名称: NRIP3 Rabbit Polyclonal Antibody
产品货号: APRab14898




研究背景 (Background)


NRIP3 (Nuclear Receptor Interacting Protein 3) is a Protein Coding gene. Gene Ontology (GO) annotations related to NRIP3 include aspartic-type endopeptidase activity. An important paralog of NRIP3 is NRIP2. proteolysis,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from COLO cells, using NRIP3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using NRIP3 Polyclonal Antibody

注意事项 (Note)

For research use only .