

产品名称: NSE2 Rabbit Polyclonal Antibody
产品货号: APRab14911



产品概述 (Summary)

产品名称 (Production Name)	NSE2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAM84B
别名 (Alternative Names)	FAM84B; BCMP101; NSE2; Protein FAM84B; Breast cancer membrane protein 101; Protein NSE2
基因 ID (Gene ID)	157638.0
蛋白 ID (SwissProt ID)	Q96KN1.The antiserum was produced against synthesized peptide derived from human FAM84B. AA range:11-60

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: NSE2 Rabbit Polyclonal Antibody
产品货号: APRab14911



研究背景 (Background)


similarity:Belongs to the FAM84 family.,tissue specificity:Expressed in esophageal squamous cell

carcinomas.,similarity:Belongs to the FAM84 family.,tissue specificity:Expressed in esophageal squamous cell carcinomas.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling; DNA / RNA; DNA Damage & Repair; Homologous Recomb; Cell Biology; Proteolysis / Ubiquitin; Proteasome / Ubiquitin; Sumo

图片 (Image Data)



Western blot analysis of lysates from Jurkat and HUVEC cells, using FAM84B Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of Jurkat cells using NSE2 Polyclonal Antibody

注意事项 (Note)

For research use only .