

产品名称: NTF2 Rabbit Polyclonal Antibody
产品货号: AP Rab14927



产品概述 (Summary)

产品名称 (Production Name)	NTF2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	NUTF2
别名 (Alternative Names)	NUTF2; NTF2; Nuclear transport factor 2; NTF-2; Placental protein 15; PP15
基因 ID (Gene ID)	10204.0
蛋白 ID (SwissProt ID)	P61970.The antiserum was produced against synthesized peptide derived from human NUTF2. AA range:67-116

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	15kDa

研究背景 (Background)


产品名称: NTF2 Rabbit Polyclonal Antibody
产品货号: APRab14927

This gene encodes a cytosolic factor that facilitates protein transport into the nucleus. The encoded protein is required for nuclear import of the small Ras-like GTPase, Ran which is involved in numerous cellular processes. This protein also interacts with the nuclear pore complex glycoprotein p62. [provided by RefSeq, Apr 2016],function:Facilitates protein transport into the nucleus. Interacts with the nucleoporin p62 and with Ran. Acts at a relatively late stage of nuclear protein import, subsequent to the initial docking of nuclear import ligand at the nuclear envelope. Could be part of a multicomponent system of cytosolic factors that assemble at the pore complex during nuclear import.,similarity:Contains 1 NTF2 domain.,subunit:Homodimer.,


研究领域 (Research Area)

Signal Transduction; Protein Trafficking; Nuclear Import / Export

图片 (Image Data)




Western blot analysis of NUTF2 Antibody. The lane on the right is blocked with the NUTF2 peptide.




Immunohistochemistry analysis of paraffin-embedded human prostate, using NUTF2 Antibody. The lane on the right is blocked with the NUTF2 peptide.


产品名称: NTF2 Rabbit Polyclonal Antibody
产品货号: APRab14927



Western blot analysis of the lysates from HT-29 cells using NTF2 antibody.



Western Blot analysis of HeLa cells using NTF2 Polyclonal Antibody



Western Blot analysis of HeLa cells using NTF2 Polyclonal Antibody

注意事项 (Note)

For research use only .