

**产品名称: NUFP2 Rabbit Polyclonal Antibody**  
**产品货号: APRab14974**



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | NUFP2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,ELISA                         |
| 种属反应性 (Reactivity)     | Human,Mouse                      |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG                                                                                      |
| 克隆 (Clonality)      | Polyclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化方式 (Purification) | Affinity purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | NUFIP2 KIAA1321 PIG1                                                         |
| 别名 (Alternative Names) |                                                                              |
| 基因 ID (Gene ID)        | 57532.0                                                                      |
| 蛋白 ID (SwissProt ID)   | Q7Z417.Synthesized peptide derived from human protein . at AA range: 180-260 |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 76kDa                                |

## 研究背景 (Background)

function: Binds RNA.,PTM: Phosphorylated upon DNA damage, probably by ATM or ATR.,subcellular location: Distribution is


**产品名称: NUFP2 Rabbit Polyclonal Antibody**  
**产品货号: APRab14974**



cell cycle-modulated, being cytoplasmic in the G2/M phase and accumulating in nucleus during the G1 phase.,subunit:Interacts with FMRP.,function:Binds RNA.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subcellular location:Distribution is cell cycle-modulated, being cytoplasmic in the G2/M phase and accumulating in nucleus during the G1 phase.,subunit:Interacts with FMRP.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western Blot analysis of HeLa lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

## 注意事项 (Note)

For research use only .