

产品名称: O13A1 Rabbit Polyclonal Antibody  
产品货号: APRab15024



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | O13A1 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | WB,ELISA                         |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                  |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG                                                                                      |
| 克隆 (Clonality)      | Polyclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化方式 (Purification) | Affinity purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                              |
|------------------------|------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | OR13A1                                                                       |
| 别名 (Alternative Names) |                                                                              |
| 基因 ID (Gene ID)        | 79290.0                                                                      |
| 蛋白 ID (SwissProt ID)   | Q8NGR1.Synthesized peptide derived from human protein . at AA range: 190-270 |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 36kDa                                |

## 研究背景 (Background)

Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception

**产品名称: O13A1 Rabbit Polyclonal Antibody**  
**产品货号: APRab15024**




of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],caution:It is uncertain whether Met-1 or Met-20 is the initiator.,function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family.,

## 研究领域 (Research Area)

Olfactory transduction;

## 图片 (Image Data)



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night

## 注意事项 (Note)

For research use only .