

**产品名称: OAS2 Rabbit Polyclonal Antibody**  
**产品货号: APRab15078**



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | OAS2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA             |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                 |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                 |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | OAS2                                                                                                                            |
| 别名 (Alternative Names) | OAS2; 2'-5'-oligoadenylate synthase 2; (2-5')oligo(A) synthase 2; 2-5A synthase 2; p69 OAS / p71 OAS; p69OAS / p71OAS           |
| 基因 ID (Gene ID)        | 4939.0                                                                                                                          |
| 蛋白 ID (SwissProt ID)   | P29728.The antiserum was produced against synthesized peptide derived from the N-terminal region of human OAS2. AA range:61-110 |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 82kDa                                                                   |

产品名称: OAS2 Rabbit Polyclonal Antibody  
产品货号: APRab15078




## 研究背景 (Background)

2'-5'-oligoadenylate synthetase 2(OAS2) Homo sapiens This gene encodes a member of the 2-5A synthetase family, essential proteins involved in the innate immune response to viral infection. The encoded protein is induced by interferons and uses adenosine triphosphate in 2'-5'-specific nucleotidyl transfer reactions to synthesize 2'-5'-oligoadenylates (2-5As). These molecules activate latent RNase L, which results in viral RNA degradation and the inhibition of viral replication. The three known members of this gene family are located in a cluster on chromosome 12. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008],catalytic activity: Binds double-stranded RNA and polymerizes ATP into PPP(A2'P5'A)N oligomers, which activate the latent RNase L that, when activated, cleaves single-stranded RNAs.,function: May play a role in mediating resistance to virus infection, control of cell growth, differentiation, and apoptosis.,induction: By interferons.,similarity: Belongs to the 2-5A synthetase family.,subcellular location: Associated with different subcellular fractions such as mitochondrial, nuclear, and rough/smooth microsomal fractions.,subunit: Homodimer.,


## 研究领域 (Research Area)

DNA / RNA; RNA Processing; Epigenetics and Nuclear Signaling; Chromatin Binding Proteins; DNA / RNA binding; Microbiology; Organism; Virus; RNA Virus; ssRNA positive strand virus; SARS Coronavirus; Immunology; Innate Immunity; Cytokines; Interferons

## 图片 (Image Data)




Western Blot analysis of HeLa cells using OAS2 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000



**产品名称: OAS2 Rabbit Polyclonal Antibody**  
**产品货号: APRab15078**



Immunohistochemical analysis of paraffin-embedded human-uterus, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-uterus, antibody was diluted at 1:100

#### **注意事项 (Note)**

For research use only .