

产品名称: Olfactory receptor 5W2 Rabbit Polyclonal Antibody
产品货号: APRab15304



产品概述 (Summary)

产品名称 (Production Name)	Olfactory receptor 5W2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OR5W2
别名 (Alternative Names)	OR5W2; OR5W2P; OR5W3P; Olfactory receptor 5W2; Olfactory receptor 5W3; Olfactory receptor OR11-155
基因 ID (Gene ID)	390148.0
蛋白 ID (SwissProt ID)	Q8NH69. The antiserum was produced against synthesized peptide derived from human OR5W2. AA range: 158-207

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

产品名称: Olfactory receptor 5W2 Rabbit Polyclonal Antibody
产品货号: AP Rab15304




研究背景 (Background)


olfactory receptor family 5 subfamily W member 2(OR5W2) Homo sapiens Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using OR5W2 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from HeLa cells, using OR5W2 Antibody. The lane on the right is blocked with the

产品名称: Olfactory receptor 5W2 Rabbit Polyclonal
Antibody
产品货号: AP Rab15304



Western blot analysis of the lysates from 293 cells using OR5W2 antibody.



Western Blot analysis of 293T-UV cells using Olfactory receptor 5W2 Polyclonal Antibody diluted at 1 : 500

注意事项 (Note)

For research use only .