

产品名称: Olfactory receptor 8J3 Rabbit Polyclonal Antibody
产品货号: APRab15332



产品概述 (Summary)

产品名称 (Production Name)	Olfactory receptor 8J3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OR8J3
别名 (Alternative Names)	OR8J3; Olfactory receptor 8J3; Olfactory receptor OR11-173
基因 ID (Gene ID)	81168.0
蛋白 ID (SwissProt ID)	Q8NGG0. The antiserum was produced against synthesized peptide derived from human OR8J3. AA range: 232-281

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)

产品名称: Olfactory receptor 8J3 Rabbit Polyclonal Antibody
产品货号: AP Rab15332




function:Odorant receptor .,similarity:Belongs to the G-protein coupled receptor 1 family.,function:Odorant receptor .,similarity:Belongs to the G-protein coupled receptor 1 family.,


研究领域 (Research Area)

Olfactory transduction;


图片 (Image Data)



Western blot analysis of lysates from HeLa, HepG2, and HUVEC cells, using OR8J3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using OR8J3 antibody.



产品名称: Olfactory receptor 8J3 Rabbit Polyclonal Antibody
产品货号: APRab15332



注意事项 (Note)

For research use only .