

产品名称: OSR2 Rabbit Polyclonal Antibody
产品货号: APRab15522



产品概述 (Summary)

产品名称 (Production Name)	OSR2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OSR2
别名 (Alternative Names)	OSR2; Protein odd-skipped-related 2
基因 ID (Gene ID)	116039.0
蛋白 ID (SwissProt ID)	Q8N2R0.Synthesized peptide derived from OSR2 . at AA range: 60-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	35kDa

研究背景 (Background)

odd-skipped related transcription factor 2(OSR2) Homo sapiens OSR2 is a mammalian homolog of the Drosophila odd-

产品名称: OSR2 Rabbit Polyclonal Antibody
产品货号: AP Rab15522




skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).[supplied by OMIM, Mar 2008],similarity:Belongs to the Odd C2H2-type zinc-finger protein family.,similarity:Contains 5 C2H2-type zinc fingers.,

研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Transcription; Domain Families; Zinc Finger; Stem Cells; Mesenchymal Stem Cells; Osteogenesis

图片 (Image Data)



Western Blot analysis of JK cells using OSR2 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .