

产品名称: p140Cap Rabbit Polyclonal Antibody
产品货号: APRab15572



产品概述 (Summary)

产品名称 (Production Name)	p140Cap Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SRCIN1
别名 (Alternative Names)	SRCIN1; KIAA1684; P140; SRC kinase signaling inhibitor 1; SNAP-25-interacting protein; SNIP; p130Cas-associated protein; p140Cap
基因 ID (Gene ID)	80725.0
蛋白 ID (SwissProt ID)	Q9C0H9.The antiserum was produced against synthesized peptide derived from human SNIP. AA range:811-860

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	112kDa

产品名称: p140Cap Rabbit Polyclonal Antibody
产品货号: APRab15572




研究背景 (Background)


function: Delays the onset of cell spreading in the early stages of cell adhesion to fibronectin. Also involved in calcium-dependent exocytosis from PC12 cells., PTM: Tyrosine-phosphorylated in response to EGF and to cell adhesion to integrin ligands., similarity: Belongs to the p140Cap family., subcellular location: Localized to the perinuclear region, lamellopodia, cortical actin and actin stress fibers but not to focal adhesions., subunit: Interacts with the N-terminal coiled-coil region of SNAP25 (By similarity). Interacts with BCAR1/p130Cas through its C-terminal domain., function: Delays the onset of cell spreading in the early stages of cell adhesion to fibronectin. Also involved in calcium-dependent exocytosis from PC12 cells., PTM: Tyrosine-phosphorylated in response to EGF and to cell adhesion to integrin ligands., similarity: Belongs to the p140Cap family., subcellular location: Localized to the perinuclear region, lamellopodia, cortical actin and actin stress fibers but not to focal adhesions., subunit: Interacts with the N-terminal coiled-coil region of SNAP25 (By similarity). Interacts with BCAR1/p130Cas through its C-terminal domain.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of various cells using p140Cap Polyclonal Antibody diluted at 1: 500



Western Blot analysis of HeLa cells using p140Cap Polyclonal Antibody diluted at 1: 500

产品名称: p140Cap Rabbit Polyclonal Antibody
产品货号: APRab15572



注意事项 (Note)

For research use only .