

产品名称: p164-RhoGEF Rabbit Polyclonal Antibody
产品货号: AP Rab15581



产品概述 (Summary)

产品名称 (Production Name)	p164-RhoGEF Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARHGEF17 ARHGEF17; KIAA0337; TEM4; Rho guanine nucleotide exchange factor 17; 164
别名 (Alternative Names)	kDa Rho-specific guanine-nucleotide exchange factor; p164-RhoGEF; p164RhoGEF; Tumor endothelial marker 4
基因 ID (Gene ID)	9828.0
蛋白 ID (SwissProt ID)	Q96PE2. The antiserum was produced against synthesized peptide derived from human ARHGEF17. AA range: 431-480

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: p164-RhoGEF Rabbit Polyclonal Antibody
产品货号: APRab15581




研究背景 (Background)

function:Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPases.,similarity:Contains 1 DH (DBL-homology) domain.,tissue specificity:Highly expressed in the heart.,function:Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPases.,similarity:Contains 1 DH (DBL-homology) domain.,tissue specificity:Highly expressed in the heart.,

研究领域 (Research Area)

Regulation of Actin Dynamics; AMPK

图片 (Image Data)



Immunofluorescence analysis of HepG2 cells, using CARD6 Antibody diluted at 1:50. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .