

产品名称: p23 Rabbit Polyclonal Antibody
产品货号: AP Rab15590



产品概述 (Summary)

产品名称 (Production Name)	p23 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PTGES3 PTGES3; P23; TEBP; Prostaglandin E synthase 3; Cytosolic prostaglandin E2 synthase; cPGES; Hsp90 co-chaperone; Progesterone receptor complex p23; Telomerase-binding protein p23
别名 (Alternative Names)	
基因 ID (Gene ID)	10728.0
蛋白 ID (SwissProt ID)	Q15185.The antiserum was produced against synthesized peptide derived from human TEBP. AA range:79-128

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

产品名称: p23 Rabbit Polyclonal Antibody
产品货号: AP Rab15590




研究背景 (Background)

This gene encodes an enzyme that converts prostaglandin endoperoxide H₂ (PGH₂) to prostaglandin E₂ (PGE₂). This protein functions as a co-chaperone with heat shock protein 90 (HSP90), localizing to response elements in DNA and disrupting transcriptional activation complexes. Alternative splicing results in multiple transcript variants. There are multiple pseudogenes of this gene on several different chromosomes. [provided by RefSeq, Feb 2016], catalytic activity: (5Z,13E)-(15S)-9- α ,11- α -epidioxy-15-hydroxyprosta-5,13-dienoate = (5Z,13E)-(15S)-11- α ,15-dihydroxy-9-oxoprost-5,13-dienoate, function: Molecular chaperone that localizes to genomic response elements in a hormone-dependent manner and disrupts receptor-mediated transcriptional activation, by promoting disassembly of transcriptional regulatory complexes, pathway: Lipid metabolism; prostaglandin biosynthesis, similarity: Belongs to the p23/wos2 family, similarity: Contains 1 CS domain, subunit: Binds to telomerase and to the progesterone receptor,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TEBP Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from NIH/3T3 cells, treated with EGF 200ng/ml 30', using TEBP Antibody. The lane on the


产品名称: p23 Rabbit Polyclonal Antibody
产品货号: APRab15590



right is blocked with the synthesized peptide.



Western Blot analysis of various cells using p23 Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of K562 cells using p23 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .