

产品名称: p53AIP1 Rabbit Polyclonal Antibody
产品货号: APRab15649



产品概述 (Summary)

产品名称 (Production Name)	p53AIP1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TP53AIP1
别名 (Alternative Names)	TP53AIP1; p53-regulated apoptosis-inducing protein 1; p53AIP1
基因 ID (Gene ID)	63970.0
蛋白 ID (SwissProt ID)	Q9HCN2. The antiserum was produced against synthesized peptide derived from human TPIP1. AA range: 75-124

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

产品名称: p53AIP1 Rabbit Polyclonal Antibody
产品货号: APRab15649




This gene is specifically expressed in the thymus, and encodes a protein that is localized to the mitochondrion. The expression of this gene is inducible by p53, and it is thought to play an important role in mediating p53-dependent apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Oct 2011],function:May play an important role in mediating p53/TP53-dependent apoptosis.,induction:By p53/TP53.,tissue specificity:Only found to be expressed in thymus.,

研究领域 (Research Area)

p53;

图片 (Image Data)



Immunofluorescence analysis of A549 cells, using TPIP1 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .