

产品名称: PAI-2 Rabbit Polyclonal Antibody
产品货号: AP Rab15701



产品概述 (Summary)

产品名称 (Production Name)	PAI-2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	SERPINB2 SERPINB2; PAI2; PLANH2; Plasminogen activator inhibitor 2; PAI-2; Monocyte
别名 (Alternative Names)	Arg-serpin; Placental plasminogen activator inhibitor; Serpin B2; Urokinase inhibitor
基因 ID (Gene ID)	5055.0
蛋白 ID (SwissProt ID)	P05120.The antiserum was produced against synthesized peptide derived from human PAI-2. AA range:258-307

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	46kDa

产品名称: PAI-2 Rabbit Polyclonal Antibody
产品货号: APRab15701




研究背景 (Background)

function:Inhibits urokinase-type plasminogen activator. The monocyte derived PAI-2 is distinct from the endothelial cell-derived PAI-1.,PTM:The signal sequence is not cleaved.,similarity:Belongs to the serpin family. Ov-serpin subfamily.,function:Inhibits urokinase-type plasminogen activator. The monocyte derived PAI-2 is distinct from the endothelial cell-derived PAI-1.,PTM:The signal sequence is not cleaved.,similarity:Belongs to the serpin family. Ov-serpin subfamily.,


研究领域 (Research Area)

Immunology; Innate Immunity; Macrophage / Inflamm; Cell Biology; Proteolysis / Ubiquitin; Protease inhibitors; Serine protease inhibitors; SERPINS

图片 (Image Data)



Immunohistochemistry analysis of PAI-2 antibody in paraffin-embedded human brain tissue.



Western blot analysis of lysate from MCF-7 cells, using PAI-2 antibody.

注意事项 (Note)

For research use only .