

产品名称: Pbx 4 Rabbit Polyclonal Antibody
产品货号: APRab15810



产品概述 (Summary)

产品名称 (Production Name)	Pbx 4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PBX4
别名 (Alternative Names)	PBX4; Pre-B-cell leukemia transcription factor 4; Homeobox protein PBX4
基因 ID (Gene ID)	80714.0
蛋白 ID (SwissProt ID)	Q9BYU1.Synthesized peptide derived from Pbx 4 . at AA range: 220-300

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

This gene encodes a member of the pre-B cell leukemia transcription factor family. These proteins are homeobox proteins


产品名称: Pbx 4 Rabbit Polyclonal Antibody
产品货号: APRab15810



that play critical roles in embryonic development and cellular differentiation both as Hox cofactors and through Hox-independent pathways. The encoded protein contains a homeobox DNA-binding domain, but specific functions of the protein have not been determined. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, May 2011],similarity:Belongs to the TALE/PBX homeobox family.,similarity:Contains 1 homeobox DNA-binding domain.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200 (4° overnight) .
2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 45min) .

注意事项 (Note)

For research use only .