

产品名称: PEA3 Rabbit Polyclonal Antibody
产品货号: AP Rab15946



产品概述 (Summary)

产品名称 (Production Name)	PEA3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ETV4 ETV4; E1AF; PEA3; ETS translocation variant 4; Adenovirus E1A enhancer-binding protein; E1A-F; Polyomavirus enhancer activator 3 homolog; Protein PEA3
别名 (Alternative Names)	
基因 ID (Gene ID)	2118.0
蛋白 ID (SwissProt ID)	P43268. The antiserum was produced against synthesized peptide derived from human ETV4. AA range: 281-330

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: PEA3 Rabbit Polyclonal Antibody
产品货号: APRab15946




研究背景 (Background)

function:Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3',similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,function:Transcriptional activator that binds to the enhancer of the adenovirus E1A gene; the core-binding sequence is 5'[AC]GGA[AT]GT-3',similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,

研究领域 (Research Area)

Microbiology; Interspecies Interaction; Host Virus Interaction; Neuroscience; Neurology process; Neurogenesis; Epigenetics and Nuclear Signaling; Transcription; Domain Families; Ets

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ETV4 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .