

产品名称: Peroxin 11 γ Rabbit Polyclonal Antibody
产品货号: AP Rab15983



产品概述 (Summary)

产品名称 (Production Name)	Peroxin 11 γ Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PEX11G PEX11G; PEX11C; Peroxisomal membrane protein 11C; Peroxin-11C;
别名 (Alternative Names)	Peroxisomal biogenesis factor 11C; Protein PEX11 homolog gamma; PEX11-gamma
基因 ID (Gene ID)	92960.0
蛋白 ID (SwissProt ID)	Q96HA9. The antiserum was produced against synthesized peptide derived from human PEX11G. AA range: 11-60

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: Peroxin 11y Rabbit Polyclonal Antibody
产品货号: AP Rab15983




研究背景 (Background)

The protein encoded by this gene is a member of the PEX11 family. This family is reported to regulate the number and size of peroxisomes in evolutionarily distant organisms. The protein encoded by this gene may induce clustering of peroxisomes. Alternative splicing results in multiple transcript variants that encode different protein isoforms. [provided by RefSeq, Jul 2012], similarity: Belongs to the peroxin-11 family.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human pancreas, using PEX11C Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .