

产品名称: PFAAP5 Rabbit Polyclonal Antibody
产品货号: APRab16009



产品概述 (Summary)

产品名称 (Production Name)	PFAAP5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	N4BP2L2
别名 (Alternative Names)	N4BP2L2; CG005; PFAAP5; NEDD4-binding protein 2-like 2; Phosphonoformate immuno-associated protein 5
基因 ID (Gene ID)	10443.0
蛋白 ID (SwissProt ID)	Q92802.The antiserum was produced against synthesized peptide derived from human N4BP2L2. AA range:358-407

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	67kDa

产品名称: PFAAP5 Rabbit Polyclonal Antibody
产品货号: AP Rab16009




研究背景 (Background)


PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from LOVO cells, using N4BP2L2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from K562 cells using N4BP2L2 antibody.



产品名称: PFAAP5 Rabbit Polyclonal Antibody
产品货号: APRab16009



Western Blot analysis of LOVO cells using PFAAP5 Polyclonal Antibody

注意事项 (Note)

For research use only .