

产品名称: PGLYRP1 Rabbit Polyclonal Antibody

产品货号: APRab16038

产品概述 (Summary)

描述(Description)	Rabbit polyclonal Antibody
宿主(Host)	Rabbit
应用(Application)	WB,IHC,ELISA
种属反应(Reactivity)	Human,Rat,Mouse
偶联物(Conjugation)	Unconjugated
修饰(Modification)	Unmodified
同种型(Isotype)	IgG
克隆(Clonality)	Polyclonal
剂型(Form)	Liquid
保存条件(Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液(Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化 (Purification)	Affinity purification

产品应用(Application)

稀释比(Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
分子量(Molecular Weight)	21kDa

抗原信息(Antigen Information)

基因名(Gene Name)	PGLYRP1
别名(Alternative Names)	PGLYRP1; PGLYRP; PGRP; TNFSF3L; SBBI68; Peptidoglycan recognition protein 1; Peptidoglycan recognition protein short; PGRP-S
基因 ID(Gene ID)	8993.0
SwissProt ID	O75594
免疫原(Immunogen)	The antiserum was produced against synthesized peptide derived from the Internal region of human PGLYRP1. AA range:131-180

研究背景 (Background)

PGLYRP1 (Peptidoglycan Recognition Protein 1) is a Protein Coding gene. Among its related pathways are TNF Signaling (sino).

GO annotations related to this gene include peptidoglycan binding and N-acetylmuramoyl-L-alanine amidase activity. An important paralog of this gene is PGLYRP4. peptidoglycan metabolic process, polysaccharide catabolic process, polysaccharide metabolic process, aminoglycan metabolic process, aminoglycan catabolic process, glycosaminoglycan catabolic process, defense response, immune response, behavior, rhythmic behavior, circadian rhythm, macromolecule catabolic process, peptidoglycan catabolic process, detection of external stimulus, detection of biotic stimulus, response to bacterium, detection of bacterium, carbohydrate catabolic process, circadian sleep/wake cycle process, glycosaminoglycan metabolic process, defense response to bacterium, circadian sleep/wake cycle, regulation of circadian sleep/wake cycle, regulation of circadian rhythm, innate immune response, regulation of circadian sleep/wake cycle, sleep, rhythmic process, circadian behavior, regulation of behavior, defense response to Gram-positive bacterium, detection of stimulus,


研究领域 (Research Area)

Cell Biology; Apoptosis; Receptors; Associated Proteins; Neuroscience; Neurology process; Circadian Rhythm; Hormones; Microbiology; Protein; Human Protein; Defensin; Immunology; Innate Immunity; Cytokines; Cancer; Cell Death

注意事项 (Note)

For research use only.

图片 (Image Data)



Western blot analysis of HepG2 823-AV 293T AD293 HeLa lysis using PGLYRP1 antibody. Antibody was diluted at 1:1000. Secondary antibody was diluted at 1:20000

EnkiLife 优势产品推荐

产品类别	产品货号	产品名称
WB 解决方案	RA10020	2h 极速 WB 即用型全流程试剂盒
	RA10021	4h 快速 WB 即用型全流程试剂盒
	RA10042	5h 畅享版 WB 全流程试剂盒
	RA10037	校准级彩色预染蛋白 Marker (8-180kDa)
	RA10038	校准级彩色预染蛋白 Marker (10-250kDa)
	RA10039	校准级高分子彩色预染蛋白 Marker (25-400kDa)
TSA 多重荧光染色试剂盒	RA10008	TSA 双标三色多重荧光染色试剂盒 (mIHC)
	RA10009	TSA 三标四色多重荧光染色试剂盒 (mIHC)
	RA10010	TSA 四标五色多重荧光染色试剂盒 (mIHC)
	RA10011	TSA 五标六色多重荧光染色试剂盒 (mIHC)
	RA10012	TSA 六标七色多重荧光染色试剂盒 (mIHC)
IHC 检测试剂盒	RA10006	HRP Anti-Mouse/Rabbit IHC Detection System
	RA10007	Polymer-HRP Anti-Mouse/Rabbit IHC Detection System
抗体标记试剂盒	RE80004p	辣根过氧化物酶(HRP)抗体标记试剂盒
	RE80002q	Sulfo-NHS-生物素标记试剂盒
	RE80007p	Cy3 荧光素标记试剂盒
	RE80011p	Fluor488 荧光素标记试剂盒
	RE80017p	Fluor750 荧光素标记试剂盒
	RE80005p	藻红蛋白(R-PE) 抗体快速标记试剂盒
	RE80040	PE-Cy7 串联染料抗体快速标记试剂盒
稳转细胞系构建服务 (免费赠送全膜 WB 验证)	TS-0001	过表达稳转细胞系构建
	TS-0002	敲低稳转细胞系构建
	TS-0003	敲除细胞系构建

限于篇幅，更多产品详见官网：<https://www.enkilife.cn>