

产品名称: PGRMC2 Rabbit Polyclonal Antibody
产品货号: AP Rab16047



产品概述 (Summary)

产品名称 (Production Name)	PGRMC2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PGRMC2
别名 (Alternative Names)	PGRMC2; DG6; PMBP; Membrane-associated progesterone receptor component 2; Progesterone membrane-binding protein; Steroid receptor protein DG6
基因 ID (Gene ID)	10424.0
蛋白 ID (SwissProt ID)	O15173.The antiserum was produced against synthesized peptide derived from human PGRMC2. AA range:174-223

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	25kDa

产品名称: PGRMC2 Rabbit Polyclonal Antibody
产品货号: APRab16047




研究背景 (Background)

domain:The cytochrome b5 heme-binding domain lacks the conserved iron-binding His residues at positions 137 and 161.,function:Receptor for steroids .,similarity:Belongs to the cytochrome b5 family. MAPR subfamily.,similarity:Contains 1 cytochrome b5 heme-binding domain.,domain:The cytochrome b5 heme-binding domain lacks the conserved iron-binding His residues at positions 137 and 161.,function:Receptor for steroids .,similarity:Belongs to the cytochrome b5 family. MAPR subfamily.,similarity:Contains 1 cytochrome b5 heme-binding domain.,

研究领域 (Research Area)

Signal Transduction; Growth Factors/Hormones; Hormones

图片 (Image Data)



Western Blot analysis of various cells using PGRMC2 Polyclonal Antibody diluted at 1 : 1000

注意事项 (Note)

For research use only .