

产品名称: PHI-1 Rabbit Polyclonal Antibody
产品货号: AP Rab16071



产品概述 (Summary)

产品名称 (Production Name)	PHI-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP1R14B
别名 (Alternative Names)	PPP1R14B; PLCB3N; PNG; Protein phosphatase 1 regulatory subunit 14B; Phospholipase C-beta-3 neighbouring gene protein
基因 ID (Gene ID)	26472.0
蛋白 ID (SwissProt ID)	Q96C90.Synthesized peptide derived from PHI-1 . at AA range: 50-130

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: PHI-1 Rabbit Polyclonal Antibody
产品货号: APRab16071



function:Inhibitor of PPP1CA. Has over 50-fold higher inhibitory activity when phosphorylated.,PTM:Phosphorylated primarily on Thr-57 by PKC (in vitro). An unknown Ser is also phosphorylated by PKC (in vitro).,similarity:Belongs to the PP1 inhibitor family.,tissue specificity:Ubiquitous. Expressed at low levels.,function:Inhibitor of PPP1CA. Has over 50-fold higher inhibitory activity when phosphorylated.,PTM:Phosphorylated primarily on Thr-57 by PKC (in vitro). An unknown Ser is also phosphorylated by PKC (in vitro).,similarity:Belongs to the PP1 inhibitor family.,tissue specificity:Ubiquitous. Expressed at low levels.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

注意事项 (Note)

For research use only .