

产品名称: Photomedin-1 Rabbit Polyclonal Antibody
产品货号: APRab16083



产品概述 (Summary)

产品名称 (Production Name)	Photomedin-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	OLFML2A
别名 (Alternative Names)	OLFML2A; Olfactomedin-like protein 2A; Photomedin-1
基因 ID (Gene ID)	169611.0
蛋白 ID (SwissProt ID)	Q68BL7.The antiserum was produced against synthesized peptide derived from human OLFML2A. AA range:311-360

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300
蛋白分子量 (Molecular Weight)	73kDa

研究背景 (Background)

产品名称: Photomedin-1 Rabbit Polyclonal Antibody
产品货号: APRab16083




PTM:May be cleaved at Lys-295 after secretion.,PTM:O-glycosylated but not N-glycosylated.,similarity:Contains 1 olfactomedin-like domain.,subunit:Homodimer. Binds to heparin and chondroitin sulfate E.,PTM:May be cleaved at Lys-295 after secretion.,PTM:O-glycosylated but not N-glycosylated.,similarity:Contains 1 olfactomedin-like domain.,subunit:Homodimer. Binds to heparin and chondroitin sulfate E.,


研究领域 (Research Area)

Signal Transduction; Cytoskeleton / ECM; Extracellular Matrix; ECM Proteins; Other ECM Proteins

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using OLFML2A Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from RAW264.7 cells, using OLFML2A Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Photomedin-1 Rabbit Polyclonal Antibody
产品货号: AP Rab16083



Western Blot analysis of RAW264.7 cells using Photomedin-1 Polyclonal Antibody

注意事项 (Note)

For research use only .