

产品名称: PITSLRE Rabbit Polyclonal Antibody
产品货号: APRab16166



产品概述 (Summary)

产品名称 (Production Name)	PITSLRE Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK11B
别名 (Alternative Names)	CDK11B; CDC2L1; CDK11; PITSLREA; PK58; Cyclin-dependent kinase 11B; Cell division cycle 2-like protein kinase 1; CLK-1; Cell division protein kinase 11B; Galactosyltransferase-associated protein kinase p58/GTA; PITSLRE serine/threonine-prot
基因 ID (Gene ID)	984.0
蛋白 ID (SwissProt ID)	P21127.The antiserum was produced against synthesized peptide derived from human CDC2L1. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
----------------------	---

产品名称: PITSLRE Rabbit Polyclonal Antibody
产品货号: APRab16166



蛋白分子量 (Molecular Weight) 92kDa


研究背景 (Background)

This gene encodes a member of the serine/threonine protein kinase family. Members of this kinase family are known to be essential for eukaryotic cell cycle control. Due to a segmental duplication, this gene shares very high sequence identity with a neighboring gene. These two genes are frequently deleted or altered in neuroblastoma. The protein kinase encoded by this gene can be cleaved by caspases and may play a role in cell apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014], caution: The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,


研究领域 (Research Area)

Epigenetics and Nuclear Signaling; Cell cycle; Kinases/Phosphatases; Cdk

图片 (Image Data)



Western blot analysis of lysates from HUVEC and HepG2 cells, using CDC2L1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using PITSLRE Polyclonal Antibody

产品名称: PITSLRE Rabbit Polyclonal Antibody
产品货号: APRab16166



注意事项 (Note)

For research use only .