

产品名称: PMCH Rabbit Polyclonal Antibody
产品货号: APRab16300



产品概述 (Summary)

产品名称 (Production Name)	PMCH Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PMCH MCH
别名 (Alternative Names)	Pro-MCH [Cleaved into: Neuropeptide-glycine-glutamic acid (NGE; Neuropeptide G-E); Neuropeptide-glutamic acid-isoleucine (NEI; Neuropeptide E-I); Melanin-concentrating hormone (MCH)]
基因 ID (Gene ID)	5367.0
蛋白 ID (SwissProt ID)	P20382.Synthetic peptide from human protein at AA range: 112-161

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: PMCH Rabbit Polyclonal Antibody
产品货号: APRab16300




研究背景 (Background)

pro-melanin concentrating hormone(PMCH) Homo sapiens This gene encodes a preproprotein that is proteolytically processed to generate multiple protein products. These products include melanin-concentrating hormone (MCH), neuropeptide-glutamic acid-isoleucine (NEI), and neuropeptide-glycine-glutamic acid (NGE). Melanin-concentrating hormone is a 19-amino acid neuropeptide that stimulates hunger and may additionally regulate energy homeostasis, reproductive function, and sleep. Pseudogenes of this gene have been identified on chromosome 5. [provided by RefSeq, Jul 2015],function:MCH may act as a neurotransmitter or neuromodulator in a broad array of neuronal functions directed toward the regulation of goal-directed behavior, such as food intake, and general arousal. May also have a role in spermatocyte differentiation.,PTM:Differentially processed in the brain and in peripheral organs producing two neuropeptides; NEI and MCH. A third peptide, NGE, may also be produced. Preferential processing in neurons by prohormone convertase 2 (PC2) generates NEI. MCH is generated in neurons of the lateral hypothalamic area by several prohormone convertases including PC1/3, PC2 and PC5/6.,similarity:Belongs to the melanin-concentrating hormone family.,tissue specificity:Predominantly expressed in lateral hypothalamus, also detected in pallidum, neocortex and cerebellum. Also found in thymus, brown adipose tissue, duodenum and testis (spermatogonia, early spermatocytes and Sertoli cells). No expression in peripheral blood. In brain exclusively mature MCH and NEI peptides are present. In peripheral tissues a large product, encompassing the NEI and MCH domains of the precursor, is found predominantly.,


研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-testis, antibody was diluted at 1:200

产品名称: PMCH Rabbit Polyclonal Antibody
产品货号: APRab16300



Immunohistochemical analysis of paraffin-embedded human-testis, antibody was diluted at 1:200

注意事项 (Note)

For research use only .