

产品名称: PML Rabbit Polyclonal Antibody
产品货号: APRab16306



产品概述 (Summary)

产品名称 (Production Name)	PML Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PML
别名 (Alternative Names)	PML; MYL; RNF71; TRIM19; Protein PML; Promyelocytic leukemia protein; RING finger protein 71; Tripartite motif-containing protein 19
基因 ID (Gene ID)	5371.0
蛋白 ID (SwissProt ID)	P29590.The antiserum was produced against synthesized peptide derived from human PML. AA range:11-60

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	98kDa

产品名称: PML Rabbit Polyclonal Antibody
产品货号: AP Rab16306




研究背景 (Background)

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. This phosphoprotein localizes to nuclear bodies where it functions as a transcription factor and tumor suppressor. Its expression is cell-cycle related and it regulates the p53 response to oncogenic signals. The gene is often involved in the translocation with the retinoic acid receptor alpha gene associated with acute promyelocytic leukemia (APL). Extensive alternative splicing of this gene results in several variations of the protein's central and C-terminal regions; all variants encode the same N-terminus. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008], alternative products: Additional isoforms seem to exist, disease: A chromosomal aberration involving PML may be a cause of acute promyelocytic leukemia (APL). Translocation t(15;17)(q21;q21) with RARA. The PML breakpoints (type A and type B) lie on either side of an alternatively spliced exon., function: Probable transcription factor. May play an important role in recruitment of ELF4 into PML nuclear bodies., PTM: Sumoylated on all three sites is required for nuclear body formation. Sumoylation on Lys-160 is a prerequisite for sumoylation on Lys-65. The PML-RARA fusion protein is not sumoylated., PTM: Ubiquitinated; mediated by SIAH1 or SIAH2 and leading to its subsequent proteasomal degradation., similarity: Contains 1 RING-type zinc finger., similarity: Contains 2 B box-type zinc fingers., subcellular location: Sumoylated forms localize to the PML nuclear bodies. The B1 box and the RING finger are also required for this nuclear localization. Isoforms lacking a nuclear localization signal are cytoplasmic., subunit: Interacts with SIRT1, TOPBP1, TRIM27 and TRIM69. Interacts with the C-terminus of ELF4. Interacts with Lassa virus Z protein and rabies virus phosphoprotein.,

研究领域 (Research Area)


Ubiquitin mediated proteolysis; Pathways in cancer; Acute myeloid leukemia;

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using PML Antibody. The picture on the right is blocked with the synthesized peptide.


产品名称: PML Rabbit Polyclonal Antibody
产品货号: APRab16306



Immunohistochemistry analysis of paraffin-embedded human tonsil tissue, using PML Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from A549 cells, using PML Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using PML Polyclonal Antibody diluted at 1: 1000

产品名称: PML Rabbit Polyclonal Antibody
产品货号: APRab16306



Western Blot analysis of A549 cells using PML Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .