

**产品名称: POTE14/22 Rabbit Polyclonal Antibody**  
**产品货号: APRab16381**



## 产品概述 (Summary)

|                        |                                      |
|------------------------|--------------------------------------|
| 产品名称 (Production Name) | POTE14/22 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody           |
| 宿主 (Host)              | Rabbit                               |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA                  |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                      |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | POTEG<br>POTEG; A26C2; POTE14; POTE ankyrin domain family member G; ANKRD26-like family C member 2; Prostate; ovary, testis-expressed protein on chromosome 14; POTE-14; POTEH; A26C3; ACTBL1; POTE22; POTE ankyrin domain family member H; ANKRD26-like |
| 别名 (Alternative Names) |                                                                                                                                                                                                                                                          |
| 基因 ID (Gene ID)        | 404785/23784                                                                                                                                                                                                                                             |
| 蛋白 ID (SwissProt ID)   | Q6S5H5/Q6S545/A6NI47.The antiserum was produced against synthesized peptide derived from human POTEG. AA range:401-450                                                                                                                                   |

## 产品应用 (Application)

|                      |                                                                         |
|----------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio) | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000 |
|----------------------|-------------------------------------------------------------------------|

产品名称: POTE14/22 Rabbit Polyclonal Antibody  
产品货号: APRab16381




蛋白分子量 (Molecular Weight) 67kDa

## 研究背景 (Background)


similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma, using A26C2/3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using POTE14/22 Polyclonal Antibody

## 注意事项 (Note)

For research use only .