

产品名称: POTE15 Rabbit Polyclonal Antibody
产品货号: APRab16383



产品概述 (Summary)

产品名称 (Production Name)	POTE15 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POTEB POTEB; A26B1; POTE15; POTE ankyrin domain family member B; ANKRD26-like
别名 (Alternative Names)	family B member 1; Prostate; ovary, testis-expressed protein on chromosome 15; POTE-15
基因 ID (Gene ID)	339010.0
蛋白 ID (SwissProt ID)	Q6S5H4.The antiserum was produced against synthesized peptide derived from human A26B1. AA range:481-530

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	60kDa

产品名称: POTE15 Rabbit Polyclonal Antibody
产品货号: AP Rab16383




研究背景 (Background)


POTEB (POTE Ankyrin Domain Family Member B) is a Protein Coding gene. An important paralog of this gene is ANKRD18A.

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using A26B1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of Mouse-kidney 3T3 lysis using POTE15 antibody. Antibody was diluted at 1:500

注意事项 (Note)

For research use only .