

产品名称: POTE-8 Rabbit Polyclonal Antibody
产品货号: AP Rab16384



产品概述 (Summary)

产品名称 (Production Name)	POTE-8 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	POTEA POTEA; A26A1; POTE8; POTE ankyrin domain family member A; ANKRD26-like
别名 (Alternative Names)	family A member 1; Prostate; ovary, testis-expressed protein on chromosome 8; POTE-8
基因 ID (Gene ID)	340441.0
蛋白 ID (SwissProt ID)	Q6S8J7. The antiserum was produced against synthesized peptide derived from human POTEA. AA range: 375-424

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: POTE-8 Rabbit Polyclonal Antibody
产品货号: APRab16384




研究背景 (Background)

similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,similarity:Belongs to the POTE family.,similarity:Contains 5 ANK repeats.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma, using POTE8 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .