

产品名称: PP2A-C α Rabbit Polyclonal Antibody
产品货号: APRab16396



产品概述 (Summary)

产品名称 (Production Name)	PP2A-C α Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP2CA
别名 (Alternative Names)	PPP2CA; Serine/threonine-protein phosphatase 2A catalytic subunit alpha isoform; PP2A-alpha; Replication protein C; RP-C
基因 ID (Gene ID)	5515.0
蛋白 ID (SwissProt ID)	P67775.The antiserum was produced against synthesized peptide derived from human PP2A-alpha. AA range:260-309

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	35kDa

产品名称: PP2A-C α Rabbit Polyclonal Antibody
产品货号: APRab16396




研究背景 (Background)

This gene encodes the phosphatase 2A catalytic subunit. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. This gene encodes an alpha isoform of the catalytic subunit. [provided by RefSeq, Jul 2008], catalytic activity: A phosphoprotein + H₂O = a protein + phosphate, cofactor: Binds 1 iron ion per subunit, cofactor: Binds 1 manganese ion per subunit, function: PP2A can modulate the activity of phosphorylase B kinase casein kinase 2, mitogen-stimulated S6 kinase, and MAP-2 kinase. Can dephosphorylate SV40 large T antigen and p53. Dephosphorylates SV40 large T antigen, preferentially on serine residues 120, 123, 677, and perhaps 679. The C subunit was most active, followed by the AC form, which was more active than the ABC form, and activity of all three forms was strongly stimulated by manganese, and to a lesser extent by magnesium. Dephosphorylation by the AC form, but not C or ABC form is inhibited by small T antigen, PTM: Phosphorylation of either threonine (by autophosphorylation-activated protein kinase) or tyrosine results in inactivation of the phosphatase. Auto-dephosphorylation has been suggested as a mechanism for reactivation, PTM: Reversibly methyl esterified on Leu-309. Carboxyl methylation may play a role in holoenzyme assembly. It varies during the cell cycle, similarity: Belongs to the PPP phosphatase family, similarity: Belongs to the PPP phosphatase family. PP-1 subfamily, subcellular location: In prometaphase cells, but not in anaphase cells, localizes at centromeres. During mitosis, also found at spindle poles, subunit: PP2A consists of a common heterodimeric core enzyme, composed of a 36 kDa catalytic subunit (subunit C) and a 65 kDa constant regulatory subunit (PR65 or subunit A), that associates with a variety of regulatory subunits. Proteins that associate with the core dimer include three families of regulatory subunits B (the R2/B/PR55/B55, R3/B''/PR72/PR130/PR59 and R5/B'/B56 families), the 48 kDa variable regulatory subunit, viral proteins, and cell signaling molecules. Interacts with NXN; the interaction is direct (By similarity). May indirectly interact with SGOL1, most probably through regulatory B56 subunits,

研究领域 (Research Area)

Oocyte meiosis; WNT; WNT-T CELL TGF-beta; Tight junction; Long-term depression;


图片 (Image Data)




产品名称: PP2A-C α Rabbit Polyclonal Antibody
产品货号: APRab16396




Immunofluorescence analysis of HepG2 cells, using PP2A-alpha Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using PP2A-alpha Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from A549 cells, using PP2A-alpha Antibody. The lane on the right is blocked with the synthesized peptide.




Western Blot analysis of various cells using PP2A-C α Polyclonal Antibody diluted at 1: 2000

产品名称: PP2A-C α Rabbit Polyclonal Antibody
产品货号: APRab16396



Western Blot analysis of A549 cells using PP2A-C α Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of various cells using Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .