

产品名称: PPP4R1L Rabbit Polyclonal Antibody
产品货号: APRab16434



产品概述 (Summary)

产品名称 (Production Name)	PPP4R1L Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PPP4R1L
别名 (Alternative Names)	PPP4R1L; C20orf192; Putative serine/threonine-protein phosphatase 4 regulatory subunit 1-like
基因 ID (Gene ID)	55370.0
蛋白 ID (SwissProt ID)	Q9P1A2.The antiserum was produced against synthesized peptide derived from human PPP4R1L. AA range:191-240

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	50kDa

产品名称: PPP4R1L Rabbit Polyclonal Antibody
产品货号: APRab16434




研究背景 (Background)


function: May be a regulatory subunit of serine/threonine-protein phosphatase 4, similarity: Contains 4 HEAT repeats, function: May be a regulatory subunit of serine/threonine-protein phosphatase 4, similarity: Contains 4 HEAT repeats,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO, HUVEC, HepG2, Jurkat, and COS cells, using PPP4R1L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using PPP4R1L antibody.

产品名称: PPP4R1L Rabbit Polyclonal Antibody
产品货号: APRab16434



Western Blot analysis of various cells using PPP4R1L Polyclonal Antibody

注意事项 (Note)

For research use only .