

产品名称: PROL4 Rabbit Polyclonal Antibody
产品货号: APRab16520



产品概述 (Summary)

产品名称 (Production Name)	PROL4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PRR4 LPRP PROL4
别名 (Alternative Names)	
基因 ID (Gene ID)	11272.0
蛋白 ID (SwissProt ID)	Q16378.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	14kDa

研究背景 (Background)

proline rich 4 (lacrimal)(PRR4) Homo sapiens This gene encodes a member of the proline-rich protein family that lacks a conserved repetitive domain. This protein may play a role in protective functions in the eye. Alternative splicing result in


产品名称: PROL4 Rabbit Polyclonal Antibody
产品货号: AP Rab16520



multiple transcript variants. Read-through transcription also exists between this gene and the upstream PRH1 (proline-rich protein HaeIII subfamily 1) gene. [provided by RefSeq, Feb 2011], tissue specificity: Abundantly expressed in lacrimal gland where it is found in the acinar cells but not in the intralobular ducts. Also found in the submandibular gland, the parotid and sublingual glands.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .