

产品名称: PSA Rabbit Polyclonal Antibody
产品货号: APRab16575



产品概述 (Summary)

产品名称 (Production Name)	PSA Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	KLK3
别名 (Alternative Names)	KLK3; APS; Prostate-specific antigen; PSA; Gamma-seminoprotein; Semin; Kallikrein-3; P-30 antigen; Semenogelase
基因 ID (Gene ID)	354.0
蛋白 ID (SwissProt ID)	P07288.The antiserum was produced against synthesized peptide derived from human Prostate-specific Antigen. AA range:81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	34kDa

产品名称: PSA Rabbit Polyclonal Antibody
产品货号: APRab16575




研究背景 (Background)

Kallikreins are a subgroup of serine proteases having diverse physiological functions. Growing evidence suggests that many kallikreins are implicated in carcinogenesis and some have potential as novel cancer and other disease biomarkers. This gene is one of the fifteen kallikrein subfamily members located in a cluster on chromosome 19. Its protein product is a protease present in seminal plasma. It is thought to function normally in the liquefaction of seminal coagulum, presumably by hydrolysis of the high molecular mass seminal vesicle protein. Serum level of this protein, called PSA in the clinical setting, is useful in the diagnosis and monitoring of prostatic carcinoma. Alternate splicing of this gene generates several transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008], catalytic activity: Preferential cleavage: -Tyr-|-Xaa-, function: Hydrolyzes semenogelin-1 thus leading to the liquefaction of the seminal coagulum, online information: Prostate-specific antigen entry, similarity: Belongs to the peptidase S1 family, similarity: Belongs to the peptidase S1 family. Kallikrein subfamily, similarity: Contains 1 peptidase S1 domain.,


研究领域 (Research Area)

Pathways in cancer; Prostate cancer;

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue, using Prostate-specific Antigen Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from mouse heart cells, using Prostate-specific Antigen Antibody. The lane on the right is

产品名称: PSA Rabbit Polyclonal Antibody
产品货号: APRab16575



blocked with the synthesized peptide.

注意事项 (Note)

For research use only .