

产品名称: PWWP2B Rabbit Polyclonal Antibody
产品货号: APRab16705



产品概述 (Summary)

产品名称 (Production Name)	PWWP2B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PWWP2B
别名 (Alternative Names)	PWWP2B; PWWP2; PWWP domain-containing protein 2B
基因 ID (Gene ID)	170394.0
蛋白 ID (SwissProt ID)	Q6NUJ5.The antiserum was produced against synthesized peptide derived from human PWWP2B. AA range:268-317

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	64kDa

研究背景 (Background)

产品名称: PWWP2B Rabbit Polyclonal Antibody
产品货号: AP Rab16705




PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translation N-terminally shortened.,similarity:Contains 1 PWWP domain.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,sequence caution:Translation N-terminally shortened.,similarity:Contains 1 PWWP domain.,


研究领域 (Research Area)

Cell Biology; Other Antibodies

图片 (Image Data)



Western blot analysis of lysates from LOVO cells, using PWWP2B Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of LOVO cells using PWWP2B Polyclonal Antibody

注意事项 (Note)

For research use only .