

产品名称: PX11A Rabbit Polyclonal Antibody
产品货号: APRab16706



产品概述 (Summary)

产品名称 (Production Name)	PX11A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PEX11A PEX11
别名 (Alternative Names)	
基因 ID (Gene ID)	8800.0
蛋白 ID (SwissProt ID)	O75192.Synthesized peptide derived from human protein . at AA range: 110-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	27kDa

研究背景 (Background)

This gene is a member of the PEX11 family, which is composed of membrane elongation factors involved in regulation of


产品名称: PX11A Rabbit Polyclonal Antibody
产品货号: APRab16706



peroxisome maintenance and proliferation. This gene product interacts with peroxisomal membrane protein 19 and may respond to outside stimuli to increase peroxisome abundance. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Oct 2012],function:May be involved in peroxisomal proliferation and may regulate peroxisomes division. May mediate binding of coatomer proteins to the peroxisomal membrane.,PTM:Seems not to be N-glycosylated.,similarity:Belongs to the peroxin-11 family.,subunit:Probably interacts with COPB2 and COPA (By similarity). Interacts with PEX19.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night

注意事项 (Note)

For research use only .