

产品名称: Pxmp2 Rabbit Polyclonal Antibody
产品货号: APRab16709



产品概述 (Summary)

产品名称 (Production Name)	Pxmp2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	PXMP2
别名 (Alternative Names)	PXMP2; PMP22; Peroxisomal membrane protein 2; 22 kDa peroxisomal membrane protein
基因 ID (Gene ID)	5827.0
蛋白 ID (SwissProt ID)	Q9NR77. The antiserum was produced against synthesized peptide derived from human PXMP2. AA range: 38-87

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: Pxmp2 Rabbit Polyclonal Antibody
产品货号: APRab16709




研究背景 (Background)

function: Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane., similarity: Belongs to the peroxisomal membrane protein PXMP2/4 family., subunit: Interacts with PEX19 and SIVA1., function: Seems to be involved in pore-forming activity and may contribute to the unspecific permeability of the peroxisomal membrane., similarity: Belongs to the peroxisomal membrane protein PXMP2/4 family., subunit: Interacts with PEX19 and SIVA1.,

研究领域 (Research Area)

Tags & Cell Markers; Subcellular Markers; Organelles; Peroxisome

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using PXMP2 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .