

产品名称: Rab 31 Rabbit Polyclonal Antibody
产品货号: APRab16751



产品概述 (Summary)

产品名称 (Production Name)	Rab 31 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB31
别名 (Alternative Names)	RAB31; RAB22B; Ras-related protein Rab-31; Ras-related protein Rab-22B
基因 ID (Gene ID)	11031.0
蛋白 ID (SwissProt ID)	Q13636.Synthesized peptide derived from Rab 31 . at AA range: 50-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	21kDa

研究背景 (Background)

Small GTP-binding proteins of the RAB family, such as RAB31, play essential roles in vesicle and granule targeting (Bao et

产品名称: Rab 31 Rabbit Polyclonal Antibody
产品货号: AP Rab16751




al., 2002 [PubMed 11784320]).[supplied by OMIM, Jul 2009],similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Highest expression in placenta and brain with lower levels in heart and lung. Not detected in liver, skeletal muscle, kidney or pancreas.,

研究领域 (Research Area)

Endocytosis;

图片 (Image Data)



Western Blot analysis of HeLa cells using Rab 31 Polyclonal Antibody

注意事项 (Note)

For research use only .