

产品名称: Rab 35 Rabbit Polyclonal Antibody
产品货号: APRab16754



产品概述 (Summary)

产品名称 (Production Name)	Rab 35 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB35
别名 (Alternative Names)	RAB35; RAB1C; RAY; Ras-related protein Rab-35; GTP-binding protein RAY; Ras-related protein Rab-1C
基因 ID (Gene ID)	11021.0
蛋白 ID (SwissProt ID)	Q15286.The antiserum was produced against synthesized peptide derived from human RAB35. AA range:99-148

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

产品名称: Rab 35 Rabbit Polyclonal Antibody
产品货号: AP Rab16754




研究背景 (Background)

function: Possesses GTPase activity., similarity: Belongs to the small GTPase superfamily. Rab family., subcellular location: Identified by mass spectrometry in melanosome fractions from stage I to stage IV., function: Possesses GTPase activity., similarity: Belongs to the small GTPase superfamily. Rab family., subcellular location: Identified by mass spectrometry in melanosome fractions from stage I to stage IV.,


研究领域 (Research Area)

Cell Biology; Cell Cycle; Cell Division; Cytokinesis; Signal Transduction; Signaling Pathway; G Protein Signaling; Small G Proteins; Ras Family; Protein Trafficking; Vesicle Transport; Regulation

图片 (Image Data)




Western blot analysis of RAB35 Antibody. The lane on the right is blocked with the RAB35 peptide.



Western blot analysis of the lysates from HT-29 cells using RAB35 antibody.

产品名称: Rab 35 Rabbit Polyclonal Antibody
产品货号: APRab16754



Western Blot analysis of various cells using Rab 35 Polyclonal Antibody

注意事项 (Note)

For research use only .