

产品名称: Rab 37 Rabbit Polyclonal Antibody
产品货号: AP Rab16755



产品概述 (Summary)

产品名称 (Production Name)	Rab 37 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB37
别名 (Alternative Names)	RAB37; Ras-related protein Rab-37
基因 ID (Gene ID)	326624.0
蛋白 ID (SwissProt ID)	Q96AX2.The antiserum was produced against synthesized peptide derived from human RAB37. AA range:170-219

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	26kDa

研究背景 (Background)

产品名称: Rab 37 Rabbit Polyclonal Antibody
产品货号: AP Rab16755




Rab proteins are low molecular mass GTPases that are critical regulators of vesicle trafficking. For additional background information on Rab proteins, see MIM 179508.[supplied by OMIM, Apr 2006],similarity:Belongs to the small GTPase superfamily. Rab family.,subcellular location:Secretory granules.,


研究领域 (Research Area)

Signal Transduction; Signaling Pathway; G Protein Signaling; Small G Proteins; Ras Family; Immunology; Immune System Diseases; Allergens; Innate Immunity; Mast Cells

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue, using RAB37 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HUVEC cells, using RAB37 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Rab 37 Rabbit Polyclonal Antibody
产品货号: AP Rab16755



Western blot analysis of the lysates from HepG2 cells using RAB37 antibody.

注意事项 (Note)

For research use only .