

产品名称: Rab 7L1 Rabbit Polyclonal Antibody
产品货号: AP Rab16763



产品概述 (Summary)

产品名称 (Production Name)	Rab 7L1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RAB7L1
别名 (Alternative Names)	RAB7L1; Ras-related protein Rab-7L1; Rab-7-like protein 1
基因 ID (Gene ID)	8934.0
蛋白 ID (SwissProt ID)	O14966.The antiserum was produced against synthesized peptide derived from human RAB7L1. AA range:115-164

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:50-1:300,ELISA 1:2000-1:20000
蛋白分子量 (Molecular Weight)	28kDa

研究背景 (Background)


产品名称: Rab 7L1 Rabbit Polyclonal Antibody
产品货号: AP Rab16763




similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Ubiquitous.,similarity:Belongs to the small GTPase superfamily. Rab family.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from HT-29 cells, using RAB7L1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using RAB7L1 antibody.



Western Blot analysis of various cells using Rab 7L1 Polyclonal Antibody

产品名称: Rab 7L1 Rabbit Polyclonal Antibody
产品货号: APRab16763



注意事项 (Note)

For research use only .