

产品名称: Rab L2A Rabbit Polyclonal Antibody
产品货号: APRab16764



产品概述 (Summary)

产品名称 (Production Name)	Rab L2A Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RABL2A
别名 (Alternative Names)	RABL2A; Rab-like protein 2A
基因 ID (Gene ID)	11159.0
蛋白 ID (SwissProt ID)	Q9UBK7.The antiserum was produced against synthesized peptide derived from human RABL2A. AA range:179-228

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	30-34kDa

研究背景 (Background)


产品名称: Rab L2A Rabbit Polyclonal Antibody
产品货号: AP Rab16764




This gene is a member of the RAB gene family which belongs to the RAS GTPase superfamily. The proteins in the family of RAS-related signaling molecules are small GTP-binding proteins that play important roles in the regulation of exocytotic and endocytotic pathways. This gene maps to the site of an ancestral telomere fusion event and may be a subtelomeric gene. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015], similarity: Belongs to the small GTPase superfamily. Rab family.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of various cells using Rab L2A Polyclonal Antibody diluted at 1: 500



Western Blot analysis of COLO205 cells using Rab L2A Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .